

BECKMAN USER DEFINED REAGENT CARTRIDGES (CX[®] SERIES 3, 4, 5, 7)

RINSING PROCEDURES

For Our LowVolume Users

Items Supplied

1 - Plastic Bag
1 - FedEx Ground Return Record

1. PREPARE

Empty any reagents left in cartridge.

2. RINSE

Briefly rinse each chamber of the cartridge under D.I. tap and shake out excess water.

3. BAG

Place rinsed cartridges and caps in bag liner in cartridge receptacle. The receptacle will hold up to 12 cartridges.

4. BOX

Once box is full, enclose a completed FedEx Ground pickup record.

5. SHIP

Seal box with preprinted return address on outside flap and call FedEx Ground at 1-800-937-4587 for pickup.

